ASPIRATOR TIPS

PRIME | ONE

Akzenta International SA, CH-6830 Chiasso, www.akzenta.ch

Non sterile!
Reprocessing and
sterilize before use!

DCWA0/AKZ
Stand 05 / 2026

GB Information on reprocessing according to DIN EN ISO 17664

With the aspirator tip, you have received a sophisticated and functional medical device.

1. Application and use

The aspirator tip fits handpieces with hose fittings according to DIN 13938 (ø 16 mm). The aspirator tip is not suitable for fitting to a non-standardised outer cone, such as e.g. Siemens Siro 1 to V or Adec. This may cause severe damage to the aspirator tip.

2. Product life:

The aspirator tip is made of high-quality plastic. It is autoclavable. The end of the product life is generally determined by the results of wear and tear. Frequent reprocessing does not affect the performance of the aspirator tip. Based on standard use, the aspirator tip can be reprocessed 100 times (standard value) in the instrument autoclave.

3. Reprocessing

3.1 General information on reprocessing

The aspirator tip has an easy-to-clean surface with no sharp edges and corners.

For cleaning/disinfection, instrument disinfection is recommended as cleaning / disinfectant agent.

Complete moistening of the aspirator tip is essential during the disinfection process.

Compliance with the current legal requirements for reprocessing medical devices has to be ensured (see e.g., www.rki.de).

Please observe the information and instructions for use of the manufacturer of the cleaning / disinfection unit and the manufacturer of the cleaning / disinfectant agent.

The aspirator tip can be sterilised in the autoclave at a maximum temperature of **134° C** without problem as long as **any metal contact is avoided**.

Please observe the following instructions for the sterilisation process:

1. Line the sterilising trays with filter paper.
2. Do not place instruments on top of, or immediately adjacent to aspirator tips.
3. The aspirator tips should not be flush with the tray rim.
4. Do not slide trays loaded with aspirator tips into tracks close to the heating rods (use middle track).
5. Avoid contact of the aspirator tips with the inner surface of the autoclave because this heats up to much higher temperatures than the sterilisation chamber.
6. Aspirator tips should not be sterilised alone. To ensure a more even temperature distribution within the autoclave, it is recommended to sterilise instruments at the same time.

Should, despite of the above mentioned measures, a deformation of the aspirator tips be observed, it is recommended to check the temperature of the autoclave.

3.2 Storage after use and transport to the processing site

Place the aspirator tips immediately after use on the patient in an instrument tray filled with a suitable cleaning / disinfectant agent (alkaline, aldehyde-free). This soaking prevents drying of residues. Reprocessing of the aspirator tips should be carried out no later than one hour after use. The transport of the aspirator tips to the processing site should be carried out in a closed instrument tray.

3.3 Cleaning and disinfection

According to the recommendations of the Robert Koch Institute (RKI), preference should be given to automated processing. Therefore, automated processing will be described first, followed by manual processing as an alternative option.

3.3.1 Automated processing

Equipment:

- Cleaning / disinfection unit
- Alkaline cleaning/disinfectant agent
- Instrument stand

Method:

- Remove aspirator tips from the instrument tray immediately before automated processing and clean thoroughly under running water so that no residues of the cleaning / disinfectant agent will be transferred into the cleaning / disinfection unit.
- Place aspirator tips in a suitable instrument stand.
- Place the instrument stand into the cleaning / disinfection unit in such manner that the spray jet is directly aimed at the aspirator tips.
- Fill cleaning/disinfectant agent into the cleaning / disinfection unit as directed by manufacturer of the agent and unit, respectively, and start the programme.
- After completion of the programme cycle, remove the aspirator tips from the cleaning / disinfection unit and dry them preferably with compressed air according to RKI recommendations. With regards to instrument stands, ensure that difficult-to-reach areas are dried as well.
- Carry out visual control regarding integrity and cleanness with suitable magnifying device. If traces of contamination can be identified on the aspirator tip after reprocessing, the cleaning / disinfection process must be repeated until no more contamination is visible.
- Final thermal disinfection of the not packed aspirator tips in the autoclave using suitable stands or sieve trays and observing the information provided under 3.1.

3.3.2 Alternative manual processing

Equipment:

- Cleaning brush
- Cleaning / disinfectant agent (alkaline, aldehyde-free)
- Instrument tray

Method:

- Remove aspirator tip from the instrument tray and thoroughly rinse clean under running water. Completely remove adherent contaminations with the cleaning brush, both internal and externally.
- Rinse aspirator tip under running water.
- Place aspirator tip in instrument tray filled with suitable cleaning / disinfectant agent. Ensure complete moistening of the aspirator tip.
- Observe manufacturer's information regarding concentration and application time.
- At the end of the application time, rinse the aspirator tip thoroughly with water.
- Dry the aspirator tip preferably with compressed air according to RKI recommendations
- Carry out visual control regarding integrity and cleanness with suitable magnifying device. If traces of contamination can be identified on the aspirator tip after reprocessing, the cleaning / disinfection process must be repeated until no more contamination is visible.
- Final thermal disinfection of the not packed aspirator tips in the autoclave using suitable stands or sieve trays and observing the information provided under 3.1.

3.4 Control and function test

After reprocessing, the integrity of the aspirator tips should be checked visually.

In case of mechanical damage, including scratches, bite marks, porosity, tears and other signs of wear and tear, the aspirator tips must be replaced by new ones.

3.5 Transport and storage

Ensure that the aspirator tips are transported and stored in a dust-free and humidity-free environment protected from recontamination.

3.6 Final information

The manufacturer ensures that the above mentioned methods are suitable for reprocessing aspirator tips.

It is the processor's responsibility to ensure that the actual processing with the equipment, materials and staff used in the reprocessing facility will achieve the desired results.

For this end, validation and routine monitoring of the procedure is usually required. Likewise, any deviation from the supplied instructions made by the processor should be carefully evaluated regarding their effectiveness and possible negative consequences.

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